

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501085

MAPK11 Mouse Monoclonal Antibody [Clone ID: OTI1C2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1C2
Applications: IF, WB

Recommend Dilution: WB 1:2000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MAPK11 (NP_002742) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 41.4 kDa

Gene Name: mitogen-activated protein kinase 11

Database Link: NP 002742 Entrez Gene 5600 Human

Background: The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as

an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation, and development. This kinase is most closely related to p38 MAP kinase, both of which can be activated by proinflammatory cytokines and environmental stress. This kinase is activated through its phosphorylation by MAP kinase kinases (MKKs), preferably by MKK6. Transcription

factor ATF2/CREB2 has been shown to be a substrate of this kinase.

Synonyms: p38-2; P38B; p38Beta; P38BETA2; PRKM11; SAPK2; SAPK2B

Protein Families: Druggable Genome, Protein Kinase







Protein Pathways:

Amyotrophic lateral sclerosis (ALS), Epithelial cell signaling in Helicobacter pylori infection, Fc epsilon RI signaling pathway, GnRH signaling pathway, Leukocyte transendothelial migration, MAPK signaling pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, Progesterone-mediated oocyte maturation, RIG-I-like receptor signaling pathway, T cell receptor signaling pathway, Toll-like receptor signaling pathway

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MAPK11 ([RC206583], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MAPK11. Positive lysates [LY400971] (100ug) and [LC400971] (20ug) can be purchased separately from OriGene.




Anti-MAPK11 mouse monoclonal antibody (TA501085) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MAPK11 ([RC206583]).





Immunofluorescent staining of HeLa cells using anti-MAPK11 mouse monoclonal antibody (TA501085).



Immunofluorescent staining of HT29 cells using anti-MAPK11 mouse monoclonal antibody (TA501085).