

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501090

BTK Mouse Monoclonal Antibody [Clone ID: OTI12D4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI12D4

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BTK (NP_000052) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 76.1 kDa

Gene Name: Bruton tyrosine kinase

Database Link: NP 000052 Entrez Gene 695 Human

Background: The protein encoded by this gene plays a crucial role in B-cell development. Mutations in this

gene cause X-linked agammaglobulinemia type 1, which is an immunodeficiency

characterized by the failure to produce mature B lymphocytes, and associated with a failure of

Ig heavy chain rearrangement.

Synonyms: AGMX1; AT; ATK; BPK; IMD1; PSCTK1; XLA


Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: B cell receptor signaling pathway, Fc epsilon RI signaling pathway, Primary immunodeficiency







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BTK ([RC211582], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BTK. Positive lysates [LY424947] (100ug) and [LC424947] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 2 different tissues and Jurkat cells by using anti-BTK monoclonal antibody (1:500).




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-BTK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501090, Dilution 1:50)






Immunohistochemical staining of paraffinembedded liver tissue within the normal limits using anti-BTK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501090, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of liver tissue using anti-BTK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501090, Dilution 1:50)



Anti-BTK mouse monoclonal antibody (TA501090) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY BTK ([RC211582]).





HEK293T cells transfected with either pCMV6-ENTRY BTK ([RC211582]) (Red) or empty vector control plasmid (Blue) were immunostained with anti-BTK mouse monoclonal (TA501090), and then analyzed by flow cytometry.