

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA501242**

## **H6PD Mouse Monoclonal Antibody [Clone ID: OTI3H5]**

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI3H5

Applications: FC, IF, WB

Recommend Dilution: WB 1:1000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human H6PD (NP\_004276) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.98 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 88.7 kDa

**Gene Name:** hexose-6-phosphate dehydrogenase/glucose 1-dehydrogenase

Database Link: NP 004276 Entrez Gene 9563 Human

**Background:** There are 2 forms of glucose-6-phosphate dehydrogenase. G form is X-linked and H form,

encoded by this gene, is autosomally linked. This H form shows activity with other hexose-6-phosphates, especially galactose-6-phosphate, whereas the G form is specific for glucose-6-phosphate. Both forms are present in most tissues, but H form is not found in red cells.


Synonyms: CORTRD1; G6PDH; GDH

**Protein Pathways:** Metabolic pathways, Pentose phosphate pathway







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY H6PD ([RC209890], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-H6PD. Positive lysates [LY401369] (100ug) and [LC401369] (20ug) can be purchased separately from OriGene.




Anti-H6PD mouse monoclonal antibody (TA501242) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY H6PD ([RC209890]).



Flow cytometric Analysis of Hela cells, using anti-H6PD antibody (TA501242), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).





Flow cytometric Analysis of Jurkat cells, using anti-H6PD antibody (TA501242), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).