

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501383

XPNPEP1 Mouse Monoclonal Antibody [Clone ID: OTI1E3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1E3

Applications: IF, IHC, WB

Recommend Dilution: WB: 1:200 - 1:1000, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human XPNPEP1 (NP_065116) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 74.6 kDa

Gene Name: X-prolyl aminopeptidase 1

Database Link: NP 065116 Entrez Gene 7511 Human

Background: This gene encodes the cytosolic form of a metalloaminopeptidase that catalyzes the cleavage

of the N-terminal amino acid adjacent to a proline residue. The gene product may play a role

in degradation and maturation of tachykinins, neuropeptides, and peptide hormones.

Alternative splicing results in multiple transcript variants.


Synonyms: APP1; SAMP; XPNPEP; XPNPEPL; XPNPEPL1

Protein Families: Druggable Genome, Protease







Product images:



HEK293T cells were transfected with pCMV6-ENTRY XPNPEP1 ([RC203535], Left lane) or the pCMV6-ENTRY control (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-XPNPEP1.





Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-XPNPEP1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).



Immunohistochemical staining of paraffinembedded Human thyroid tissue within the normal limits using anti-XPNPEP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501383, Dilution 1:50)







Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-XPNPEP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501383, Dilution 1:50)

Anti-XPNPEP1 mouse monoclonal antibody (TA501383) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY XPNPEP1 ([RC203535]).