

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501429

Coronin 1a (CORO1A) Mouse Monoclonal Antibody [Clone ID: OTI1A5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1A5

Applications: FC, IF, WB

Recommend Dilution: WB 1:500~2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CORO1A(NP_009005) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 50.8 kDa **Gene Name:** coronin 1A

Database Link: NP 009005 Entrez Gene 11151 Human

Background: This gene encodes a member of the WD repeat protein family. WD repeats are minimally

conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Alternative splicing results in multiple transcript variants. A related pseudogene has been defined on chromosome 16.


[provided by RefSeq]

Synonyms: CLABP; CLIPINA; HCORO1; IMD8; p57; TACO







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CORO1A ([RC210753], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CORO1A. Positive lysates [LY402085] (100ug) and [LC402085] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-CORO1A monoclonal antibody.




Anti-CORO1A mouse monoclonal antibody (TA501429) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY CORO1A ([RC210753]).






HEK293T cells transfected with either [RC210753] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-CORO1A antibody (TA501429), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-CORO1A antibody (TA501429), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-CORO1A antibody (TA501429), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).