

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501616

PBP (PEBP1) Mouse Monoclonal Antibody [Clone ID: OTI2D4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2D4

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PEBP1 (NP_002558) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.35 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 21.0 kDa

Gene Name: phosphatidylethanolamine binding protein 1

Database Link: NP 002558 Entrez Gene 5037 Human

Background: Binds ATP, opioids and phosphatidylethanolamine. Has lower affinity for phosphatidylinositol


and phosphatidylcholine. Serine protease inhibitor which inhibits thrombin, neuropsin and

chymotrypsin but not trypsin,tissue type plasminogen activator and elastase


Synonyms: HCNP; HCNPpp; HEL-210; HEL-S-34; HEL-S-96; PBP; PEBP-1; RKIP



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PEBP1 ([RC206355], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PEBP1. Positive lysates [LY400913] (100ug) and [LC400913] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-PEBP1 monoclonal antibody.