

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA501645**

## HSPC150 (UBE2T) Mouse Monoclonal Antibody [Clone ID: OTI2F5]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2F5
Applications: IF, WB

Recommend Dilution: WB 1:500, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human UBE2T (NP\_054895) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.33 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 22.3 kDa

Gene Name: ubiquitin conjugating enzyme E2 T

Database Link: NP 054895 Entrez Gene 29089 Human

**Background:** The covalent conjugation of ubiquitin to proteins regulates diverse cellular pathways and

proteins. Ubiquitin is transferred to a target protein through a concerted action of a ubiquitin-activating enzyme (E1), a ubiquitin-conjugating enzyme (E2), such as UBE2T, and a

ubiquitin ligase (E3) (Machida et al., 2006 [PubMed 16916645]). [supplied by OMIM]


**Synonyms:** FANCT; HSPC150; PIG50

**Protein Families:** Druggable Genome







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBE2T ([RC200748], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBE2T. Positive lysates [LY402285] (100ug) and [LC402285] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-UBE2T monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).



Anti-UBE2T mouse monoclonal antibody (TA501645) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY UBE2T ([RC200748]).