

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501720

KIST (UHMK1) Mouse Monoclonal Antibody [Clone ID: OTI2H4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2H4

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human UHMK1 (NP_787062) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.82 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 46.4 kDa

Gene Name: U2AF homology motif kinase 1

Database Link: NP 787062 Entrez Gene 127933 Human

Background: The gene encodes a serine/threonine protein kinase that promotes cell cycle progression

through G1 by phosphorylation of the cyclin-dependent kinase inhibitor 1B (p27Kip1), which causes nuclear export and degradation. The encoded protein is also thought to function in the adult nervous system and the gene has been associated with schizophrenia. Alternative

splicing results in multiple transcript variants. [provided by RefSeq]


Synonyms: KIS; KIST; P-CIP2

Protein Families: Druggable Genome, Protein Kinase







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UHMK1 ([RC214962], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UHMK1. Positive lysates [LY406187] (100ug) and [LC406187] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-UHMK1 monoclonal antibody.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-UHMK1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501720)






Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-UHMK1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501720)



Anti-UHMK1 mouse monoclonal antibody (TA501720) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY UHMK1 ([RC214962]).



Flow cytometric Analysis of Hela cells, using anti-UHMK1 antibody (TA501720), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).