

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501778

Doublecortin (DCX) Mouse Monoclonal Antibody [Clone ID: OTI4A3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4A3

Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DCX (NP_835365) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.73 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 40.4 kDa

Gene Name: doublecortin

Database Link: NP 835365 Entrez Gene 1641 Human

Background: This gene encodes a member of the doublecortin family. The protein encoded by this gene is

a cytoplasmic protein and contains two doublecortin domains, which bind microtubules. In the developing cortex, cortical neurons must migrate over long distances to reach the site of their final differentiation. The encoded protein appears to direct neuronal migration by regulating the organization and stability of microtubules. In addition, the encoded protein interacts with LIS1, the regulatory gamma subunit of platelet activating factor acetylhydrolase, and this interaction is important to proper microtubule function in the developing cortex. Mutations in this gene cause abnormal migration of neurons during development and disrupt

the layering of the cortex, leading to epilepsy, mental retardation, subcortical band heterotopia ('double cortex' syndrome) in females and lissencephaly ('smooth brain' syndrome) in males. Multiple transcript variants encoding different isoforms have been

found for this gene. [provided by RefSeq]







Synonyms: DBCN; DC; LISX; SCLH; XLIS

Protein Families: Druggable Genome


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DCX ([RC221891], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DCX. Positive lysates [LY403598] (100ug) and [LC403598] (20ug) can be purchased separately from OriGene.




Immunofluorescent staining of 293T cells transiently transfected by pCMV6-ENTRY DCX ([RC221891]) using anti-DCX mouse monoclonal antibody (TA501778, green, upper left; merged, upper right). Cell nuclei were stained with DAPI (blue, lower left). 293T cells transfected with empty vector plasmid served as negative control (merged, lower right) (1:100).




Flow cytometric Analysis of Jurkat cells, using anti-DCX antibody (TA501778), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).





Western blot analysis of extracts (35ug) from 2 tissues by using anti-DCX monoclonal antibody (1:250).



Living HEK293T cells transfected with either [RC221891] plasmid (red) or empty vector (blue) were immunostained by anti-DCX antibody (TA501778) or isotype control antibody (green), and then analyzed by flow cytometry (1:100).