

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA501810**

## Adenylosuccinate Lyase (ADSL) Mouse Monoclonal Antibody [Clone ID: OTI2D10]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2D10

Applications: FC, IHC, WB

**Recommend Dilution:** WB 1:500~2000, IHC 1:150, FLOW 1:100

Reactivity: Human, Monkey, Rat, Dog

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ADSL (NP\_000017) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.55 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 54.7 kDa

**Gene Name:** adenylosuccinate lyase

Database Link: NP 000017 Entrez Gene 315150 RatEntrez Gene 474499 DogEntrez Gene 709259

MonkeyEntrez Gene 158 Human

**Background:** Adenylsuccinate lyase is involved in both de novo synthesis of purines and formation of

adenosine monophosphate from inosine monophosphate. It catalyzes two reactions in AMP biosynthesis: the removal of a fumarate from succinylaminoimidazole carboxamide (SAICA) ribotide to give aminoimidazole carboxamide ribotide (AICA) and removal of fumarate from adenylosuccinate to give AMP. Adenylosuccinase deficiency results in succinylpurinemic autism, psychomotor retardation, and , in some cases, growth retardation associated with muscle wasting and epilepsy. Two transcript variants encoding different isoforms have been

found for this gene. [provided by RefSeq]

Synonyms: AMPS; ASASE; ASL

Protein Families: Druggable Genome







**Protein Pathways:** 

Alanine, aspartate and glutamate metabolism, Metabolic pathways, Purine metabolism


## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ADSL ([RC200524], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ADSL. Positive lysates [LY424970] (100ug) and [LC424970] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ADSL monoclonal antibody.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)






Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)






Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)



Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-ADSL mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501810)



HEK293T cells transfected with either [RC200524] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ADSL antibody (TA501810), and then analyzed by flow cytometry.