

# OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA501863**

### AK3 Mouse Monoclonal Antibody [Clone ID: OTI6G4]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI6G4

**Applications:** FC, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human AK3 (NP\_057336) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.39 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 25.4 kDa

**Gene Name:** adenylate kinase 3

Database Link: NP 057366 Entrez Gene 50808 Human

Synonyms: AK3L1; AK6; AKL3L; AKL3L1; FIX


**Protein Families:** Druggable Genome

**Protein Pathways:** Pyrimidine metabolism







#### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AK3 ([RC204408], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AK3. Positive lysates [LY414079] (100ug) and [LC414079] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-AK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501863)



HEK293T cells transfected with either [RC204408] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-AK3 antibody (TA501863), and then analyzed by flow cytometry.

