

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501942

OTUB2 Mouse Monoclonal Antibody [Clone ID: OTI6F7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6F7

Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human OTUB2 (NP_075601) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.42 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 27.0 kDa

Gene Name: OTU deubiquitinase, ubiquitin aldehyde binding 2

Database Link: NP 075601 Entrez Gene 78990 Human

Background: Otubains are deubiquitylating cysteine proteases (DUBs; see MIM 602519) that belong to the

ovarian tumor (OTU) protein superfamily. Like other DUBs, otubains cleave proteins precisely

at the ubiquitin (UB; see MIM 191339)-protein bond. [supplied by OMIM]


Synonyms: C14orf137; OTB2; OTU2

Protein Families: Protease







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY OTUB2 ([RC209650], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-OTUB2. Positive lysates [LY402958] (100ug) and [LC402958] (20ug) can be purchased separately from OriGene.



Anti-OTUB2 mouse monoclonal antibody (TA501942) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY OTUB2 ([RC209650]).



HEK293T cells transfected with either [RC209650] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-OTUB2 antibody (TA501942), and then analyzed by flow cytometry.