

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA502037**

## PANK3 Mouse Monoclonal Antibody [Clone ID: OTI3C4]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI3C4
Applications: FC, WB

**Recommend Dilution:** WB 1:500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PANK3 (NP\_078870) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 40.9 kDa

**Gene Name:** pantothenate kinase 3

Database Link: NP 078870 Entrez Gene 79646 Human

**Background:** This gene encodes a protein belonging to the pantothenate kinase family. Pantothenate

kinase is a key regulatory enzyme in the biosynthesis of coenzyme A (CoA) in bacteria and mammalian cells. It catalyzes the first committed step in the universal biosynthetic pathway leading to CoA and is itself subject to regulation through feedback inhibition by CoA. This

family member is expressed most abundantly in the liver. [provided by RefSeq].

COMPLETENESS: complete on the 3' end.

**Synonyms:** DKFZp686O22247; FLJ12899; MGC16863


**Protein Families:** Druggable Genome

**Protein Pathways:** Metabolic pathways, Pantothenate and CoA biosynthesis







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PANK3 ([RC203696], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PANK3. Positive lysates [LY411210] (100ug) and [LC411210] (20ug) can be purchased separately from OriGene.




HEK293T cells transfected with either [RC203696] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PANK3 antibody (TA502037), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-PANK3 antibody (TA502037), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).





Flow cytometric Analysis of Jurkat cells, using anti-PANK3 antibody (TA502037), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).