

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502288

Butyrylcholinesterase (BCHE) Mouse Monoclonal Antibody [Clone ID: OTI4C12]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4C12

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BCHE (NP_000046) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 68.2 kDa

Gene Name: butyrylcholinesterase

Database Link: NP 000046 Entrez Gene 590 Human

Background: Mutant alleles at the BCHE locus are responsible for suxamethonium sensitivity. Homozygous

persons sustain prolonged apnea after administration of the muscle relaxant

suxamethonium in connection with surgical anesthesia. The activity of pseudocholinesterase in the serum is low and its substrate behavior is atypical. In the absence of the relaxant, the

homozygote is at no known disadvantage. [provided by RefSeq]


Synonyms: CHE1; CHE2; E1

Protein Families: Druggable Genome, Transmembrane





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BCHE ([RC202198], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BCHE. Positive lysates [LY424950] (100ug) and [LC424950] (20ug) can be purchased separately from OriGene.