

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA502509**

## NORE1 (RASSF5) Mouse Monoclonal Antibody [Clone ID: OTI1H2]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1H2
Applications: FC, WB

**Recommend Dilution:** WB 1:500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RASSF5 (NP\_872606) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.41 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 30.2 kDa

Gene Name: Ras association domain family member 5

Database Link: NP 872606 Entrez Gene 83593 Human

**Background:** This gene is a member of the Ras association domain family. It functions as a tumor

suppressor, and is inactivated in a variety of cancers. The encoded protein localizes to centrosomes and microtubules, and associates with the GTP-activated forms of Ras, Rap1, and several other Ras-like small GTPases. The protein regulates lymphocyte adhesion and suppresses cell growth in response to activated Rap1 or Ras. Multiple transcript variants

encoding different isoforms have been found for this gene. [provided by RefSeq]

Synonyms: Maxp1; NORE1; NORE1A; NORE1B; RAPL; RASSF3


**Protein Families:** Druggable Genome

**Protein Pathways:** Leukocyte transendothelial migration, Non-small cell lung cancer, Pathways in cancer







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RASSF5 ([RC203854], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RASSF5. Positive lysates [LY405440] (100ug) and [LC405440] (20ug) can be purchased separately from OriGene.




HEK293T cells transfected with either [RC203854] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-RASSF5 antibody (TA502509), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-RASSF5 antibody (TA502509), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).





Flow cytometric Analysis of Jurkat cells, using anti-RASSF5 antibody (TA502509), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).