

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502537

HAGHL Mouse Monoclonal Antibody [Clone ID: OTI3H2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3H2
Applications: FC, WB

Recommend Dilution: WB 1:500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HAGHL (NP_115680) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.18 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 31.1 kDa

Gene Name: hydroxyacylglutathione hydrolase like

Database Link: NP 115680 Entrez Gene 84264 Human


Synonyms: MGC2605

Protein Pathways: Pyruvate metabolism






Product images:





HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HAGHL ([RC200832], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HAGHL. Positive lysates [LY410246] (100ug) and [LC410246] (20ug) can be purchased separately from OriGene.




HEK293T cells transfected with either [RC200832] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HAGHL antibody (TA502537), and then analyzed by flow cytometry.









Flow cytometric Analysis of Jurkat cells, using anti-HAGHL antibody (TA502537), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).