

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA502602**

## **GBA3** Mouse Monoclonal Antibody [Clone ID: OTI1F1]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI1F1

**Applications:** FC, IHC, WB

**Recommend Dilution:** WB 1:500~2000, IHC 1:150, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 1-150 and 370-469 of

human GBA3 (NP\_066024) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1.2 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 53.5 kDa

**Gene Name:** glucosylceramidase beta 3 (gene/pseudogene)

Database Link: NP 066024 Entrez Gene 57733 Human

Background: GBA3, or cytosolic beta-glucosidase (EC 3.2.1.21), is a predominantly liver enzyme that

efficiently hydrolyzes beta-D-glucoside and beta-D-galactoside, but not any known physiologic beta-glycoside, suggesting that it may be involved in detoxification of plant glycosides (de Graaf et al., 2001 [PubMed 11389701]). GBA3 also has significant neutral

glycosylceramidase activity (EC 3.2.1.62), suggesting that it may be involved in a

nonlysosomal catabolic pathway of glucosylceramide metabolism (Hayashi et al., 2007

[PubMed 17595169]). [supplied by OMIM]


Synonyms: CBG; CBGL1; GLUC; KLRP

**Protein Pathways:** Cyanoamino acid metabolism, Starch and sucrose metabolism







## **Product images:**



Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)






Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)






Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)






Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-GBA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502602)





HEK293T cells transfected with either [RC211035] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-GBA3 antibody (TA502602), and then analyzed by flow cytometry.




Flow cytometric Analysis of Jurkat cells, using anti-GBA3 antibody (TA502602), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).










HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GBA3 ([RC211035], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GBA3. Positive lysates [LY402815] (100ug) and [LC402815] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-GBA3 monoclonal antibody.