

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502671

HAND1 Mouse Monoclonal Antibody [Clone ID: OTI1G10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1G10
Applications: FC, IF, WB

Recommend Dilution: WB 1:500, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HAND1 (NP_004812) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.82 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 23.4 kDa

Gene Name: heart and neural crest derivatives expressed 1

Database Link: NP 004812 Entrez Gene 9421 Human

Background: The protein encoded by this gene belongs to the basic helix-loop-helix family of transcription

factors. This gene product is one of two closely related family members, the HAND proteins, which are asymmetrically expressed in the developing ventricular chambers and play an essential role in cardiac morphogenesis. Working in a complementary fashion, they function in the formation of the right ventricle and aortic arch arteries, implicating them as mediators of congenital heart disease. In addition, it has been suggested that this transcription factor

may be required for early trophoblast differentiation. [provided by RefSeq]


Synonyms: bHLHa27; eHand; Hxt; Thing1

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HAND1 ([RC200513], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HAND1. Positive lysates [LY417727] (100ug) and [LC417727] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 9 Human tissue by using anti-HAND1 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon).




Anti-HAND1 mouse monoclonal antibody (TA502671) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HAND1 ([RC200513]).





Flow cytometric Analysis of Hela cells, using anti-HAND1 antibody (TA502671), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-HAND1 antibody (TA502671), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).