

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502679

RNF156 (MGRN1) Mouse Monoclonal Antibody [Clone ID: OTI1A11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1A11
Applications: FC, WB

Recommend Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MGRN1 (NP_056061) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.52 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 63 kDa

Gene Name: mahogunin ring finger 1

Database Link: NP 056061 Entrez Gene 23295 Human

Background: Mahogunin (MGRN1) is a C3HC4 RING-containing protein with E3 ubiquitin ligase activity in

vitro. [supplied by OMIM]


Synonyms: RNF156

Protein Pathways: Ubiquitin mediated proteolysis







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MGRN1 ([RC208284], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MGRN1. Positive lysates [LY414714] (100ug) and [LC414714] (20ug) can be purchased separately from OriGene.




HEK293T cells transfected with either [RC208284] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-MGRN1 antibody (TA502679), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-MGRN1 antibody (TA502679), (Red), compared to a nonspecific negative control antibody, (Blue).





Flow cytometric Analysis of Jurkat cells, using anti-MGRN1 antibody (TA502679), (Red), compared to a nonspecific negative control antibody, (Blue).