

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA502690**

## RNF156 (MGRN1) Mouse Monoclonal Antibody [Clone ID: OTI2A7]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2A7

**Applications:** FC, IF, WB

**Recommend Dilution:** WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MGRN1 (NP\_056061) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.53 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 63 kDa

**Gene Name:** mahogunin ring finger 1

Database Link: NP 056061 Entrez Gene 23295 Human

**Background:** Mahogunin (MGRN1) is a C3HC4 RING-containing protein with E3 ubiquitin ligase activity in

vitro. [supplied by OMIM]


Synonyms: RNF156

**Protein Pathways:** Ubiquitin mediated proteolysis







## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MGRN1 ([RC208284], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MGRN1. Positive lysates [LY414714] (100ug) and [LC414714] (20ug) can be purchased separately from OriGene.





Anti-MGRN1 mouse monoclonal antibody (TA502690) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MGRN1 ([RC208284]).



HEK293T cells transfected with either [RC208284] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-MGRN1 antibody (TA502690), and then analyzed by flow cytometry.







Flow cytometric Analysis of Jurkat cells, using anti-MGRN1 antibody (TA502690), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).