

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502929

S6K1 (RPS6KB1) Mouse Monoclonal Antibody [Clone ID: OTI6B2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6B2

Applications: IF, IHC, WB

Recommend Dilution: WB 1:500~2000, IHC 1:150, IF 1:100

Reactivity: Human, Monkey, Mouse, Dog

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RPS6KB1(NP_003152) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.7 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 59 kDa

Gene Name: ribosomal protein S6 kinase B1

Database Link: NP 003152 Entrez Gene 72508 MouseEntrez Gene 480580 DogEntrez Gene 711979

MonkeyEntrez Gene 6198 Human

Background: This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine

kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates several residues of the S6 ribosomal protein. The kinase activity of this protein leads to an increase in protein synthesis and cell proliferation. Amplification of the region of DNA encoding this gene and overexpression of this kinase are seen in some breast cancer cell lines. Alternate translational start sites have been described and alternate transcriptional splice variants have been observed but have not been thoroughly characterized. [provided by

RefSeq]

Synonyms: p70 S6KA; p70(S6K)-alpha; p70-alpha; p70-S6K; PS6K; S6K; S6K-beta-1; S6K1; STK14A

Protein Families: Druggable Genome, Protein Kinase







Protein Pathways:

Acute myeloid leukemia, ErbB signaling pathway, Fc gamma R-mediated phagocytosis, Insulin signaling pathway, mTOR signaling pathway, TGF-beta signaling pathway


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RPS6KB1 ([RC217324], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RPS6KB1. Positive lysates [LY401097] (100ug) and [LC401097] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-RPS6KB1 monoclonal antibody.




Western blot analysis of extracts (10ug) from 9 Human tissue by using anti-RPS6KB1 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon).





Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-RPS6KB1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502929)



Anti-RPS6KB1 mouse monoclonal antibody (TA502929) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY RPS6KB1 ([RC217324]).