

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503001

OSGEP Mouse Monoclonal Antibody [Clone ID: OTI8B3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8B3

Applications: IF, IHC, WB

Recommend Dilution: WB 1:500, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-335 of human OSGEP

(NP_060277) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.58 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 36.2 kDa

Gene Name: O-sialoglycoprotein endopeptidase


Database Link: NP 060277 Entrez Gene 55644 Human

Synonyms: GCPL1; KAE1; OSGEP1; PRSMG1


Protein Families: Druggable Genome, Protease



Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY OSGEP ([RC207089], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-OSGEP. Positive lysates [LY402616] (100ug) and [LC402616] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-OSGEP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503001)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-OSGEP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503001)




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-OSGEP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503001)




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-OSGEP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503001)






Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-OSGEP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503001)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-OSGEP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503001)



Anti-OSGEP mouse monoclonal antibody (TA503001) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY OSGEP ([RC207089]).





Immunofluorescent staining of HeLa cells using anti-OSGEP mouse monoclonal antibody (TA503001).