

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA503215**

## **GPR48 (LGR4) Mouse Monoclonal Antibody [Clone ID: OTI5B8]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5B8
Applications: FC, WB

**Recommend Dilution:** WB 1:500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 230-540 of human

LGR4 (NP\_060960) produced in HEK293T Cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.78 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Gene Name:** leucine rich repeat containing G protein-coupled receptor 4

Database Link: NP 060960 Entrez Gene 55366 Human

**Background:** G protein-coupled receptors (GPCRs) play key roles in a variety of physiologic functions.

Members of the leucine-rich GPCR (LGR) family, such as GPR48, have multiple N-terminal leucine-rich repeats (LRRs) and a 7-transmembrane domain (Weng et al., 2008 [PubMed


18424556]). [supplied by OMIM]

Synonyms: BNMD17; GPR48


**Protein Families:** Druggable Genome, GPCR, Transmembrane



### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LGR4 ([RC221345], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LGR4. Positive lysates [LY413017] (100ug) and [LC413017] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC221345] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-LGR4 antibody (TA503215), and then analyzed by flow cytometry.