

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503275

NUDT9 Mouse Monoclonal Antibody [Clone ID: OTI5B5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5B5

Applications: WB

Recommend Dilution: WB 1:500~2000

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NUDT9(NP_076952) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.73 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 38.9 kDa

Gene Name: nudix hydrolase 9

Database Link: NP 076952 Entrez Gene 53343 Human


Synonyms: NUDT10

Protein Families: Druggable Genome, Ion Channels: Other

Protein Pathways: Purine metabolism



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NUDT9 ([RC222894], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NUDT9. Positive lysates [LY402975] (100ug) and [LC402975] (20ug) can be purchased separately from OriGene.