

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503304

XPNPEP3 Mouse Monoclonal Antibody [Clone ID: OTI5G2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5G2

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human, Monkey, Mouse, Rat, Dog

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human XPNPEP3(NP_071381) produced in HEK293

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.7 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 56.9 kDa

Gene Name: X-prolyl aminopeptidase 3

Database Link: NP 071381 Entrez Gene 321003 MouseEntrez Gene 685823 RatEntrez Gene 481237

<u>DogEntrez Gene 705778 MonkeyEntrez Gene 63929 Human</u>

Background: The protein encoded by this gene belongs to the family of X-pro-aminopeptidases that utilize

a metal cofactor, and remove the N-terminal amino acid from peptides with a proline residue in the penultimate position. This protein has been shown to localize to the mitochondria of renal cells, and have a role in ciliary function. Mutations in this gene are associated with nephronophthisis-like nephropathy-1. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene, however, expression of some of these

isoforms in vivo is not known.


Synonyms: APP3; ICP55; NPHPL1

Protein Families: Druggable Genome, Protease







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY XPNPEP3 ([RC200888], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-XPNPEP3. Positive lysates [LY402903] (100ug) and [LC402903] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-XPNPEP3 monoclonal antibody.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-XPNPEP3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503304)






Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-XPNPEP3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503304)



Anti-XPNPEP3 mouse monoclonal antibody (TA503304) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY XPNPEP3 ([RC200888]).



HEK293T cells transfected with either [RC200888] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-XPNPEP3 antibody (TA503304), and then analyzed by flow cytometry.