

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503316

LGR5 Mouse Monoclonal Antibody [Clone ID: OTI2A2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A2

Applications: FC, IF, IHC, IP, WB

Recommend Dilution: WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100, IP 2ug/500ul

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 250-550 of human

LGR5 (NP_003658) produced in HEK293T Cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 99.8 kDa

Database Link: NP 003658 Entrez Gene 14160 Mouse

Synonyms: FEX; GPR49; GPR67; GRP49; HG38

Protein Families: Druggable Genome, GPCR, Transmembrane



Product images:




Figure from citation: Positive correlation between LSD1, c-Myc, β -catenin and LGR5 expression in human colorectal tumor tissues. Expression levels of LSD1, LGR5, β -catenin and c-Myc in consecutive sections from normal colon, adenoma and CRC tumor tissues. Dilution: 1:100 <u>View Citation</u>





Figure from citation: Immunohistochemical analysis of Lgr5 expression in human normal colon biopsies; Lgr5+ cells are localized at the crypt base, where bona fide stem-like cells home; magnification is indicated in the boxes. Dilution: 1:200 <u>View Citation</u>










Figure from citation: Immunohistochemical staining of Lgr5 in normal gastric tissues. The normal gastric tissue, 40×magnified (A) and 400×magnified (B). Yellowspots represent the Lgr5 positive cells that are localized only at base of the glands in normal gastric mucosa. The arrows in the pictures show the immunoreactive cells stained in yellow color. View Citation




Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-LGR5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503316)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-LGR5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503316)






Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-LGR5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503316)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-LGR5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503316)




Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-LGR5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503316)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-LGR5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503316)




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LGR5 ([RC212825], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LGR5. Positive lysates [LY401213] (100ug) and [LC401213] (20ug) can be purchased separately from OriGene.


Immunoprecipitation (IP) of LGR5 by using TrueMab monoclonal anti-LGR5 antibodies (Negative control: IP without adding anti-LGR5 antibody.). For each experiment, 500ul of DDK tagged LGR5 overexpression lysates (at 1:5 dilution with HEK293T lysate), 2ug of anti-LGR5 antibody and 20ul (0.1mg) of goat anti-mouse conjugated magnetic beads were mixed and incubated overnight. After extensive wash to remove any non-specific binding, the immunoprecipitated products were analyzed with rabbit anti-DDK polyclonal antibody.





Anti-LGR5 mouse monoclonal antibody (TA503316) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY LGR5 ([RC212825]).










Figure from citation: Immunocytochemistry of epithelium stem cell phenotype. Isolated dental epithelium following 5-day culture, overwhelmingly expressed Lgr5. Dilution: 1:100 <u>View Citation</u>



Immunofluorescent staining of HT29 cells using anti-LGR5 mouse monoclonal antibody (TA503316, green). Actin filaments were labeled with TRITC-phalloidin (red), and nuclear with DAPI (blue).





HEK293T cells transfected with either [RC212825] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-LGR5 antibody (TA503316), and then analyzed by flow cytometry.




Figure from citation: Cytofluorimetric analysis of Lgr5 expression by PKH high cells over 4 weeks of culture. Data are expressed as percent mean values (6 SD) of 5 consecutive experiments. Dilution: 1:400 <u>View Citation</u>




Figure from citation: Flow cytometry of the SB cells. The G4 region was further divided into P1 and P2. Lgr5, a stem cell marker, was expressed by 32% of the P2 population. Black: staining with Lgr5; red: staining with the isotype control. <u>View Citation</u>






Figure from citation: [Nature] Cells expressing the stem cell marker Lgr5 was detected by immunostaining of hCECs in 2D culture. <u>View Citation</u>

Pubmed id: 27398792 Memo: Human; IF