

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503343

COASY Mouse Monoclonal Antibody [Clone ID: OTI2H1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2H1

Applications: FC, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human COASY(NP_079509) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.59 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 62.1 kDa

Gene Name: Coenzyme A synthase

Database Link: NP 079509 Entrez Gene 80347 Human

Biosynthesis of coenzyme A (CoA) from pantothenic acid (vitamin B5) is an essential universal

pathway in prokaryotes and eukaryotes. COASY is a bifunctional enzyme that catalyzes the 2

last steps in CoA synthesis. These activities are performed by 2 separate enzymes, phosphopantetheine adenylyltransferase (PPAT; EC 2.7.7.3) and dephospho-CoA kinase (DPCK; EC 2.7.1.24), in prokaryotes (Daugherty et al., 2002 [PubMed 11923312]). [supplied by

OMIM]


Synonyms: DPCK; NBIA6; NBP; pOV-2; PPAT; UKR1

Protein Pathways: Metabolic pathways, Pantothenate and CoA biosynthesis







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY COASY ([RC220733], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-COASY. Positive lysates [LY403068] (100ug) and [LC403068] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-COASY mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503343)



Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-COASY mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503343)





HEK293T cells transfected with either [RC220733] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-COASY antibody (TA503343), and then analyzed by flow cytometry.