

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA503372**

### TFIIF (GTF2F1) Mouse Monoclonal Antibody [Clone ID: OTI1D5]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1D5

Applications: FC, IF, WB

Recommend Dilution: WB 1:1000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human GTF2F1(NP\_002087) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.88 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 58.1 kDa

Gene Name: general transcription factor IIF subunit 1

Database Link: NP 002087 Entrez Gene 2962 Human


Synonyms: BTF4; RAP74; TF2F1; TFIIF

**Protein Families:** Druggable Genome, Transcription Factors


**Protein Pathways:** Basal transcription factors



## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GTF2F1 ([RC201294], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GTF2F1. Positive lysates [LY400767] (100ug) and [LC400767] (20ug) can be purchased separately from OriGene.




Anti-GTF2F1 mouse monoclonal antibody (TA503372) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY GTF2F1 ([RC201294]).





HEK293T cells transfected with either [RC201294] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-GTF2F1 antibody (TA503372), and then analyzed by flow cytometry.



Flow cytometric Analysis of Jurkat cells, using anti-GTF2F1 antibody (TA503372), (Red), compared to a nonspecific negative control antibody, (Blue).