

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

Product datasheet for TA503423

Lipoamide Dehydrogenase (DLD) Mouse Monoclonal Antibody [Clone ID: OTI4D5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4D5

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DLD(NP_000099) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.72 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 50.1 kDa

Gene Name: dihydrolipoamide dehydrogenase

Database Link: NP 000099 Entrez Gene 1738 Human

Background: This gene encodes the L protein of the mitochondrial glycine cleavage system. The L protein,

also named dihydrolipoamide dehydrogenase, is also a component of the pyruvate dehydrogenase complex, the alpha-ketoglutarate dehydrogenase complex, and the

branched-chain alpha-keto acide dehydrogenase complex. Mutations in this gene have been

identified in patients with E3-deficient maple syrup urine disease and lipoamide

dehydrogenase deficiency. [provided by RefSeq]

Synonyms: DLDD; DLDH; E3; GCSL; LAD; PHE3

Protein Families: Druggable Genome

Protein Pathways: Citrate cycle (TCA cycle), Glycine, serine and threonine metabolism, Glycolysis /


Gluconeogenesis, Metabolic pathways, Pyruvate metabolism, Valine, leucine and isoleucine


degradation






Product images:








HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DLD ([RC200639], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DLD. Positive lysates [LY400041] (100ug) and [LC400041] (20ug) can be purchased separately from OriGene.


Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and DLD-Knockout HeLa cells (KO, Cat# [LC832718]) were separated by SDS-PAGE and immunoblotted with anti-DLD monoclonal antibody TA503423 (1:500). Then the blotted membrane was stripped and reprobed with anti-PCNA antibody as a loading control.

Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-DLD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503423)






Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-DLD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503423)



Anti-DLD mouse monoclonal antibody (TA503423) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY DLD ([RC200639]).



HEK293T cells transfected with either [RC200639] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DLD antibody (TA503423), and then analyzed by flow cytometry.