

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503449

PMVK Mouse Monoclonal Antibody [Clone ID: OTI2F9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F9

Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PMVK(NP_006547) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.32 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 21.8 kDa

Gene Name: phosphomevalonate kinase

Database Link: NP 006547 Entrez Gene 10654 Human

Background: This gene encodes a peroxisomal enzyme that catalyzes the conversion of mevalonate 5-

phosphate into mevalonate 5-diphosphate, the fifth reaction of the cholesterol biosynthetic pathway. Studies in rat show that the message level and the enzyme activity of this protein is regulated by sterol, and that this regulation is coordinated with 3-hydroxy-3-methylglutaryl coenzyme A reductase, the rate-limiting enzyme of cholesterol biosynthesis. [provided by

RefSeq, Sep 2011]

Synonyms: HUMPMKI; PMK; PMKA; PMKASE; POROK1


Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Terpenoid backbone biosynthesis







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PMVK ([RC202867], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PMVK. Positive lysates [LY416571] (100ug) and [LC416571] (20ug) can be purchased separately from OriGene.



Anti-PMVK mouse monoclonal antibody (TA503449) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PMVK ([RC202867]).



HEK293T cells transfected with either [RC202867] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PMVK antibody (TA503449), and then analyzed by flow cytometry.