

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503687

ZNF397 Mouse Monoclonal Antibody [Clone ID: OTI1F3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1F3
Applications: IF, WB

Recommend Dilution: WB 1:2000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ZNF397(NP_115723) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 30.9 kDa

Gene Name: zinc finger protein 397

Database Link: NP 115723 Entrez Gene 84307 Human

Background: This gene encodes a protein with a N-terminal SCAN domain, and the longer isoform contains

nine C2H2-type zinc finger repeats in the C-terminal domain. The protein localizes to centromeres during interphase and early prophase, and different isoforms can repress or activate transcription in transfection studies. Multiple transcript variants encoding different isoforms have been found for this gene. Additional variants have been described, but their

biological validity has not been determined. [provided by RefSeq]


Synonyms: ZNF47; ZSCAN15

Protein Families: Transcription Factors






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF397 ([RC203025], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF397. Positive lysates [LY410184] (100ug) and [LC410184] (20ug) can be purchased separately from OriGene.



Anti-ZNF397 mouse monoclonal antibody (TA503687) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ZNF397 ([RC203025]).