

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA503689

# ERAS Mouse Monoclonal Antibody [Clone ID: OTI2C6]

## **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTI2C6
Applications:	IF, WB
<b>Recommend Dilution:</b>	WB 1:2000, IF 1:100
Reactivity:	Human
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ERAS(NP_853510) produced in HEK293T cell.
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Predicted Protein Size:	25.1 kDa
Gene Name:	ES cell expressed Ras
Database Link:	<u>NP_853510 Entrez Gene 3266 Human</u>
Synonyms:	HRAS2; HRASP
Protein Families:	Druggable Genome




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ERAS ([RC210965], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ERAS. Positive lysates [LY405694] (100ug) and [LC405694] (20ug) can be purchased separately from OriGene.



Anti-ERAS mouse monoclonal antibody (TA503689) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ERAS ([RC210965]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US