

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA503713**

## **DAND5 Mouse Monoclonal Antibody [Clone ID: OTI3A3]**

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI3A3
Applications: FC, IF

Recommend Dilution: IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 23-189 of human

DAND5(NP\_689867) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 18 kDa

**Gene Name:** DAN domain BMP antagonist family member 5

Database Link: NP 689867 Entrez Gene 199699 Human

Background: This gene encodes a member of the BMP (bone morphogenic protein) antagonist family. Like

BMPs, BMP antagonists contain cystine knots and typically form homo- and heterodimers. The CAN (cerberus and dan) subfamily of BMP antagonists, to which this gene belongs, is characterized by a C-terminal cystine knot with an eight-membered ring. The antagonistic effect of the secreted protein encoded by this gene is likely due to its direct binding to BMP proteins. As an antagonist of BMP, this gene may play a role in regulating organogenesis, body patterning, and tissue differentiation. In mouse, this protein has been shown to bind Nodal and to inhibit the Nodal signaling pathway which patterns left/right body asymmetry.


[provided by RefSeq]

Synonyms: CER2; CERL2; CKTSF1B3; COCO; CRL2; DANTE; GREM3; SP1






## **Product images:**



Anti-DAND5 mouse monoclonal antibody (TA503713) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY DAND5 ([RC221703]).



HEK293T cells transfected with either [RC221703] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DAND5 antibody (TA503713), and then analyzed by flow cytometry.