

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503737

EDG2 (LPAR1) Mouse Monoclonal Antibody [Clone ID: OTI1G6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1G6
Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human LPAR1(NP_001392) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.8 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 40.9 kDa

Gene Name: lysophosphatidic acid receptor 1

Database Link: NP 001392 Entrez Gene 1902 Human

Background: The integral membrane protein encoded by this gene is a lysophosphatidic acid (LPA)

receptor from a group known as EDG receptors. These receptors are members of the G protein-coupled receptor superfamily. Utilized by LPA for cell signaling, EDG receptors mediate diverse biologic functions, including proliferation, platelet aggregation, smooth muscle contraction, inhibition of neuroblastoma cell differentiation, chemotaxis, and tumor cell invasion. Two transcript variants encoding the same protein have been identified for this

gene [provided by RefSeq, Jul

Synonyms: edg-2; EDG2; Gpcr26; GPR26; LPA1; Mrec1.3; rec.1.3; vzg-1; VZG1

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Gap junction, Neuroactive ligand-receptor interaction

