

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503772

EFEMP1 Mouse Monoclonal Antibody [Clone ID: OTI1D9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1D9
Applications: FC, WB

Recommend Dilution: WB 1:200~500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human EFEMP1(NP_004096) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 52.7 kDa

Gene Name: EGF containing fibulin extracellular matrix protein 1

Database Link: NP 004096 Entrez Gene 2202 Human

Background: This gene encodes a member of the fibulin family of extracellular matrix glycoproteins. Like

all members of this family, the encoded protein contains tandemly repeated epidermal growth factor-like repeats followed by a C-terminus fibulin-type domain. This gene is upregulated in malignant gliomas and may play a role in the aggressive nature of these tumors. Mutations in this gene are associated with Doyne honeycomb retinal dystrophy. Alternatively spliced transcript variants that encode the same protein have been described.

[provided by RefSeq]

Synonyms: DHRD; DRAD; EGF-containing fibulin-like extracellular matrix protein 1; FBLN3; FBNL; fibrillin-


like; fibulin 3; FLJ35535; MGC111353; MLVT; MTLV; S1-5


Protein Families: Secreted Protein





Product images:





HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY EFEMP1 ([RC214953], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EFEMP1. Positive lysates [LY418212] (100ug) and [LC418212] (20ug) can be purchased separately from OriGene.

Flow cytometric Analysis of Jurkat cells, using anti-EFEMP1 antibody (TA503772), (Red), compared to a nonspecific negative control antibody, (Blue).