

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503815

WDR4 Mouse Monoclonal Antibody [Clone ID: OTI4G1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4G1
Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human WDR4(NP_387510) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1.2 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 45.3 kDa

Gene Name: WD repeat domain 4

Database Link: NP 387510 Entrez Gene 10785 Human

Background: This gene encodes a member of the WD repeat protein family. WD repeats are minimally

conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene is excluded as a candidate for a form of nonsyndromic deafness (DFNB10), but is still a candidate for other disorders mapped to 21q22.3 as well as for the development of Down syndrome phenotypes. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq,


May 2012]

Synonyms: TRM82; TRMT82







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY WDR4 ([RC217569], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-WDR4. Positive lysates [LY403257] (100ug) and [LC403257] (20ug) can be purchased separately from OriGene.



Anti-WDR4 mouse monoclonal antibody (TA503815) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY WDR4 ([RC217569]).



HEK293T cells transfected with either [RC217569] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-WDR4 antibody (TA503815), and then analyzed by flow cytometry.