

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503817

HARS2 Mouse Monoclonal Antibody [Clone ID: OTI5A9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5A9

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:500, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HARS2(NP_036340) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.76 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 53.3 kDa

Gene Name: histidyl-tRNA synthetase 2, mitochondrial

Database Link: NP 036340 Entrez Gene 23438 Human

Background: Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate

amino acids. The protein encoded by this gene is an enzyme belonging to the class II family of

aminoacyl-tRNA synthetases. Functioning in the synthesis of histidyl-transfer RNA, the

enzyme plays an accessory role in the regulation of protein biosynthesis. The gene is located in a head-to-head orientation with HARS on chromosome five, where the homologous genes

share a bidirectional promoter. [provided by RefSeq]. COMPLETENESS: full length.


Synonyms: HARSL; HARSR; HO3; PRLTS2

Protein Pathways: Aminoacyl-tRNA biosynthesis







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HARS2 ([RC204925], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HARS2. Positive lysates [LY415917] (100ug) and [LC415917] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-HARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503817)




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-HARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503817)





Anti-HARS2 mouse monoclonal antibody (TA503817) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HARS2 ([RC204925]).



HEK293T cells transfected with either [RC204925] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HARS2 antibody (TA503817), and then analyzed by flow cytometry.