

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503905

ARL11 Mouse Monoclonal Antibody [Clone ID: OTI1E1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1E1

Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ARL11(NP_612459) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.7 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 21.2 kDa

Gene Name: ADP ribosylation factor like GTPase 11

Database Link: NP 612459 Entrez Gene 115761 Human

Background: This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of

proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers. [provided by


RefSeq, May 2010]

Synonyms: ARLTS1







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ARL11 ([RC203868], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ARL11. Positive lysates [LY408600] (100ug) and [LC408600] (20ug) can be purchased separately from OriGene.



Anti-ARL11 mouse monoclonal antibody (TA503905) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ARL11 ([RC203868]).



HEK293T cells transfected with either [RC203868] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ARL11 antibody (TA503905), and then analyzed by flow cytometry.