

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA504151

ARL5B Mouse Monoclonal Antibody [Clone ID: OTI2D7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2D7

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ARL5B(NP_848930) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.62 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 20.2 kDa

Gene Name: ADP ribosylation factor like GTPase 5B

Database Link: NP 848930 Entrez Gene 221079 Human

Background: ARL5B (ARL8) belongs to a family of proteins that are structurally similar to ADP-ribosylation

factors (ARFs; see MIM 103180). ARLs and ARFs are part of the RAS superfamily of regulatory


GTPases. [supplied by OMIM, Nov 2010]

Synonyms: ARL8





Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ARL5B ([RC204805], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ARL5B. Positive lysates [LY405801] (100ug) and [LC405801] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-ARL5B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504151)





Anti-ARL5B mouse monoclonal antibody (TA504151) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ARL5B ([RC204805]).



Flow cytometric Analysis of Jurkat cells, using anti-ARL5B antibody (TA504151), (Red), compared to a nonspecific negative control antibody, (Blue).