

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA504224**

## Hsp75 (TRAP1) Mouse Monoclonal Antibody [Clone ID: OTI1A1]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI1A1

Applications: IF, IHC, WB

**Recommend Dilution:** WB 1:500~2000, IHC 1:150, IF 1:100

Reactivity: Human, Monkey, Mouse, Rat, Dog

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full lengt human recombinant protein of human TRAP1(NP\_057376) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1.1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 79.9 kDa

**Gene Name:** TNF receptor associated protein 1

Database Link: NP 057376 Entrez Gene 68015 MouseEntrez Gene 287069 RatEntrez Gene 102153194

<u>DogEntrez Gene 706234 MonkeyEntrez Gene 10131 Human</u>

**Background:** HSP90 proteins are highly conserved molecular chaperones that have key roles in signal

transduction, protein folding, protein degradation, and morphologic evolution. HSP90 proteins normally associate with other cochaperones and play important roles in folding newly synthesized proteins or stabilizing and refolding denatured proteins after stress. TRAP1 is a mitochondrial HSP90 protein. Other HSP90 proteins are found in cytosol (see HSP90AA1; MIM 140571) and endoplasmic reticulum (HSP90B1; MIM 191175) (Chen et al.,

2005 [PubMed 16269234]). [supplied by OMIM]


Synonyms: HSP 75; HSP75; HSP90L; TRAP-1

**Protein Families:** Druggable Genome







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRAP1 ([RC203439], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRAP1. Positive lysates [LY402533] (100ug) and [LC402533] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-TRAP1 monoclonal antibody.




Western blot analysis of extracts (10ug) from 9 Human tissue by using anti-TRAP1 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon).






Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-TRAP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504224)




Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-TRAP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504224)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-TRAP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504224)





Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-TRAP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504224)



Anti-TRAP1 mouse monoclonal antibody (TA504224) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TRAP1 ([RC203439]).