

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA504560**

## Cystatin SA (CST2) Mouse Monoclonal Antibody [Clone ID: OTI2D6]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2D6
Applications: FC, WB

Recommend Dilution: WB 1:200~1000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human CST2(NP\_001313) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.59 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 14.3 kDa **Gene Name:** cystatin SA

Database Link: NP 001313 Entrez Gene 1470 Human

**Background:** The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences.

Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens. The type 2 cystatin proteins are a class of cysteine proteinase inhibitors found in a variety of human fluids and secretions, where they appear to provide protective functions. The cystatin

locus on chromosome 20 contains the majority of the type 2 cystatin genes and

pseudogenes. This gene is located in the cystatin locus and encodes a secreted thiol protease

inhibitor found at high levels in saliva, tears and seminal plasma. [provided by RefSeq]


Synonyms: MGC71924

**Protein Families:** Secreted Protein







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CST2 ([RC209350], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CST2. Positive lysates [LY420010] (100ug) and [LC420010] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC209350] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-CST2 antibody (TA504560), and then analyzed by flow cytometry.



Flow cytometric Analysis of Jurkat cells, using anti-CST2 antibody (TA504560), (Red), compared to a nonspecific negative control antibody, (Blue).