

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA504561

RMC1 Mouse Monoclonal Antibody [Clone ID: OTI4E4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4E4
Applications: FC, WB

Recommend Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human C18orf8(NP_037458) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.87 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 74.8 kDa

Gene Name: regulator of MON1-CCZ1

Database Link: NP 037458 Entrez Gene 29919 Human
Synonyms: C18orf8; HsT2591; Mic-1; MIC1; WDR98