

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA504588**

## TAF2 Mouse Monoclonal Antibody [Clone ID: OTI2A5]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2A5

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

**Host:** Mouse

**Isotype:** lgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 782-1121 of human

TAF2(NP\_003175) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 136.8 kDa

**Gene Name:** TATA-box binding protein associated factor 2

Database Link: NP 003175 Entrez Gene 6873 Human

**Background:** Initiation of transcription by RNA polymerase II requires the activities of more than 70

polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold

for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs. TAFs may participate in basal transcription, serve as coactivators, function in promoter recognition or modify general transcription factors (GTFs) to facilitate complex assembly and transcription initiation. This gene encodes one of the larger subunits of TFIID that is stably associated with the TFIID complex. It contributes to interactions at and downstream of the transcription initiation site, interactions that help determine transcription complex response to activators.

[provided by RefSeq]




### TAF2 Mouse Monoclonal Antibody [Clone ID: OTI2A5] - TA504588

**Synonyms:** CIF150; MRT40; TAF2B; TAFII150

**Protein Families:** Protease, Transcription Factors

**Protein Pathways:** Basal transcription factors

## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TAF2 ([RC213159], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TAF2. Positive lysates [LY418826] (100ug) and [LC418826] (20ug) can be purchased separately from OriGene.