

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA504703

RUNX1T1 Mouse Monoclonal Antibody [Clone ID: OTI1H1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1H1
Applications: IF, WB

Recommend Dilution: WB 1:2000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RUNX1T1(NP_783553) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 63 kDa

Gene Name: RUNX1 partner transcriptional co-repressor 1

Database Link: NP 783553 Entrez Gene 862 Human

Background: This gene encodes a member of the myeloid translocation gene family which interact with

DNA-bound transcription factors and recruit a range of corepressors to facilitate

transcriptional repression. The t(8;21)(q22;q22) translocation is one of the most frequent karyotypic abnormalities in acute myeloid leukemia. The translocation produces a chimeric gene made up of the 5'-region of the runt-related transcription factor 1 gene fused to the 3'-

region of this gene. The chimeric protein is thought to associate with the nuclear

corepressor/histone deacetylase complex to block hematopoietic differentiation. Alternative

splicing results in multiple transcript variants. [provided by RefSeq]

Synonyms: AML1T1; CBFA2T1; CDR; ETO; MTG8; ZMYND2


Protein Families: Transcription Factors

Protein Pathways: Acute myeloid leukemia, Pathways in cancer







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RUNX1T1 ([RC222426], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RUNX1T1. Positive lysates [LY403575] (100ug) and [LC403575] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 7 different cell lines by using anti-RUNX1T1 monoclonal antibody at 1:200 dilution.



Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-RUNX1T1 monoclonal antibody (1:200).





Anti-RUNX1T1 mouse monoclonal antibody (TA504703) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY RUNX1T1 ([RC222426]).