

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA504748**

## **DBT Mouse Monoclonal Antibody [Clone ID: OTI1G2]**

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI1G2
Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DBT(NP\_001909) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 46.4 kDa

**Gene Name:** dihydrolipoamide branched chain transacylase E2

Database Link: NP 001909 Entrez Gene 1629 Human

Background: The branched-chain alpha-keto acid dehydrogenase complex (BCKD) is an inner-

mitochondrial enzyme complex involved in the breakdown of the branched-chain amino acids isoleucine, leucine, and valine. The BCKD complex is thought to be composed of a core of 24 transacylase (E2) subunits, and associated decarboxylase (E1), dehydrogenase (E3), and regulatory subunits. This gene encodes the transacylase (E2) subunit. Mutations in this gene result in maple syrup urine disease, type 2. Alternatively spliced transcript variants have been

described, but their biological validity has not been determined. [provided by RefSeq]

Synonyms: BCATE2; BCKAD-E2; BCCADC-E2; E2; E2B


**Protein Families:** Druggable Genome

**Protein Pathways:** Metabolic pathways, Valine, leucine and isoleucine degradation







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DBT ([RC201998], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DBT. Positive lysates [LY419653] (100ug) and [LC419653] (20ug) can be purchased separately from OriGene.



Anti-DBT mouse monoclonal antibody (TA504748) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY DBT ([RC201998]).



HEK293T cells transfected with either [RC201998] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DBT antibody (TA504748), and then analyzed by flow cytometry.