

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA505094

RAB20 Mouse Monoclonal Antibody [Clone ID: OTI1F1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1F1
Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RAB20(NP_060287) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 26.1 kDa

Gene Name: RAB20, member RAS oncogene family

Database Link: NP 060287 Entrez Gene 55647 Human


Synonyms: FLJ20429

Protein Families: Druggable Genome





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RAB20 ([RC206834], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RAB20. Positive lysates [LY413526] (100ug) and [LC413526] (20ug) can be purchased separately from OriGene.