

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA505497

NKIRAS1 Mouse Monoclonal Antibody [Clone ID: OTI8B11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8B11

Applications: IF, IHC, WB

Recommend Dilution: WB 1:1000, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NKIRAS1(NP_065078) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 21.5 kDa

Gene Name: NFKB inhibitor interacting Ras like 1


Database Link: NP 065078 Entrez Gene 28512 Human

Synonyms: kappaB-Ras1; KBRAS1






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NKIRAS1 ([RC204248], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NKIRAS1. Positive lysates [LY412541] (100ug) and [LC412541] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-NKIRAS1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA505497)





Anti-NKIRAS1 mouse monoclonal antibody (TA505497) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY NKIRAS1 ([RC204248]).