

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA505684

ADAMTS8 Mouse Monoclonal Antibody [Clone ID: OTI2H5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2H5
Applications: IF, WB

Recommend Dilution: WB 1:1000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ADAMTS8(NP_008968) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 96.3 kDa

Gene Name: ADAM metallopeptidase with thrombospondin type 1 motif 8

Database Link: NP 008968 Entrez Gene 11095 Human

Background: This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with

thrombospondin motifs) protein family. Members of the family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The enzyme encoded by this gene contains two C-terminal TS motifs, and disrupts angiogenesis in vivo. A number of disorders have been mapped in the vicinity of this gene, most notably lung

neoplasms. [provided by RefSeq, Jul 2008]


Synonyms: ADAM-TS8; METH2

Protein Families: Druggable Genome, Protease, Secreted Protein






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ADAMTS8 ([RC212446], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ADAMTS8. Positive lysates [LY402079] (100ug) and [LC402079] (20ug) can be purchased separately from OriGene.



Anti-ADAMTS8 mouse monoclonal antibody (TA505684) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ADAMTS8 ([RC212446]).



Immunofluorescent staining of HeLa cells using anti-ADAMTS8 mouse monoclonal antibody (TA505684).