

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA505824**

## VTI1A Mouse Monoclonal Antibody [Clone ID: OTI1F4]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI1F4
Applications: IHC, WB

Recommend Dilution: WB 1:200~4000, IHC 1:150

Reactivity: Human, Monkey, Mouse, Rat, Dog

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human VTI1A(NP\_660207) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 25 kDa

**Gene Name:** vesicle transport through interaction with t-SNAREs 1A

Database Link: NP 660207 Entrez Gene 53611 MouseEntrez Gene 65277 RatEntrez Gene 100685447

DogEntrez Gene 696654 MonkeyEntrez Gene 143187 Human

**Synonyms:** MMDS3; MVti1; Vti1-rp2; VTI1RP2


**Protein Families:** Transmembrane

**Protein Pathways:** SNARE interactions in vesicular transport







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY VTI1A ([RC213934], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-VTI1A. Positive lysates [LY407945] (100ug) and [LC407945] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-VTI1A monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-VTI1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA505824)






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-VTI1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA505824)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-VTI1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA505824)



Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-VTI1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA505824)





Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-VTI1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA505824)