

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506094

Chromogranin A (CHGA) Mouse Monoclonal Antibody [Clone ID: OTI2F2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F2

Applications: IF, IHC, WB

Recommend Dilution: WB 1:4000, IHC 1:500, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CHGA(NP_001266) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 48.9 kDa

Gene Name: chromogranin A

Database Link: NP 001266 Entrez Gene 1113 Human

Background: The protein encoded by this gene is a member of the chromogranin/secretogranin family of

neuroendocrine secretory proteins. It is found in secretory vesicles of neurons and endocrine

cells. This gene product is a precursor to three biologically active peptides; vasostatin, pancreastatin, and parastatin. These peptides act as autocrine or paracrine negative modulators of the neuroendocrine system. Other peptides, including chromostatin, betagranin, WE-14 and GE-25, are also derived from the full-length protein. However, biological

activities for these molecules have not been shown. [provided by RefSeq, Jul 2008]


Synonyms: CGA

Protein Families: Druggable Genome, Secreted Protein







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CHGA ([RC200492], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CHGA. Positive lysates [LY400511] (100ug) and [LC400511] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-CHGA mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506094)



Immunofluorescent staining of HeLa cells using anti-CHGA mouse monoclonal antibody (TA506094).